

Ramsbottom Carbon Residue Bath

ASTM D524; IP 14; BS 2000 Part 14; BS 4451 (obs); ISO 4262; FTM 791 5002; NF T60-117



- 5 Test Wells
- Temperature range ambient to 550°C
- Digital temperature control
- Bath and Control bulb temperature displays

The Ramsbottom Carbon Residue Bath is a solid block bath with five test wells, digital temperature control, and an independent digital thermometer.

The cast iron block bath is heated by two independent heaters, a 300W disk heater in the base, and a 900W wrap-around band type heater. The controls and display for the digital temperature controller and the separate digital thermometer and thermocouple port (for measuring control bulb temperature) are mounted on the front.

SPECIFICATIONS

Part Number	10900-4
Block Temperature Range	Ambient to 550°C
Power	1300W
Voltage	110/120V, 50/60Hz 220/240V, 50/60Hz
Size (HxWxD)	40 x 42 x 42cm
Weight	33kg