



SETA REID VAPOUR PRESSURE BATH

QUALITY DESIGN, COMPACT & LOW COST

Traditional vapour pressure testing in accordance with ASTM D323 & equivalent test methods...

- FUEL BLENDING
- QUALITY CONTROL
- FUEL PERFORMANCE
- TRANSPORT & STORAGE

KEY FEATURES

- REDUCED POWER REQUIREMENTS
- GREATER EFFICIENCY – LOWER RUNNING COSTS
- INTEGRAL TEMPERATURE CONTROL
- EXCELLENT TEMPERATURE STABILITY
- REDUCED BATH HEIGHT FOR EASY INSTALLATION & REMOVAL OF VESSELS
- NO REQUIREMENT FOR EXTERNAL AIR PUMP
- SMALL FOOTPRINT
- INTEGRAL STIRRER
- 4 VESSEL SUPPORT BRACKETS INCLUDED WITH BATH



Seta RVP Bath is a compact unit that reduces floorspace requirement; it is energy efficient and easy to operate. The bath includes a thermally insulated stainless steel water tank which provides temperature precision of +/- 0.1C within the range of ambient to 75C. A digital LED controller with integral circulator ensure uniform bath temperature. An integral vessel rack allows vessels to be stored on the outside of the bath when they are not in use.



STANHOPE-SETA

Tel: +44 (0)1932 564391 Fax: +44 (0)1932 568363

E-mail: sales@stanhope-seta.co.uk Website: www.stanhope-seta.co.uk



TECHNICAL EXCELLENCE FROM STANHOPE-SETA

TEST METHODS:

ASTM D323; D1267; D4953 PROCEDURE A; IP 69; IP 410; BS2000 PART 69, 410; ISO 3007; ISO 4256; EN 12; FTMS 791B METHOD 1201; NF 12; NF M07-007; NF M41-010 (LPG); JIS K2258



Seta RVP Bath Technical Specification

Seta Product ref	22210-3
Bath	Floor standing water bath
Temperature range	ambient to 75C +/- 0.1C
Control	Integral temperature control with thermostir
Capacity	80 litres (water)
Test vessels	6 (4 if ullage chambers are fitted)
Vessel housing	Vessel rack fitted to outside of bath
Vessel supports	4 are supplied with the bath
Water Supply	Mains water inlet connection and easy access drain tap for tank emptying.
Power	2 KW
Footprint(HxWxD)	85 x 59 x 42cm
Voltage	220/240 volt 50/60 Hz

A range of Seta RVP test vessels and accessories are available:



Assembly for Gasoline



Assembly for LPG



Pressure Gauges



Pressure Calibrator



STANHOPE-SETA

London Street, Chertsey
Surrey KT16 8AP, England
Tel: +44 (0)1932 564391
Fax: +44 (0)1932 568363

E-mail: sales@stanhope-seta.co.uk
Website: www.stanhope-seta.co.uk



Distributed by

