

Seta Cloud and Pour Point Bath 11010-2

ASTM D97; ASTM D2500; ASTM D5853; IP 15; IP 219; IP 441;
BS 2000 Parts 15 & 219; ISO 3015; ISO 3016

The Cloud and Pour Point test bath determines the lowest temperature at which petroleum oils will flow and the temperature at which wax crystals form. This identifies the minimum safe operating temperature.

The bath accommodates four test jackets and a bath thermometer and comes with a stainless steel cover and a drain tap.



Specifications:

Temperature Range	9 to -69 °C
Test Positions	4
Size (HxWxD):	17 x 32 x 18 cm
Weight:	4.3 kg