

# Seta Oxi Solid Block Bath 4 Way

ASTM D525, ASTM D873, BS 2000 Parts 40 & 138, EN ISO 7536, IP 138, IP 40

- 4 position
- Dry block
- Over temperature cut out
- Small footprint

The Seta Oxi Solid Block Bath is an aluminium block bath with digital temperature control that has been designed to accept up to four oxidation vessels for carrying out tests equivalent to the IP 40 test method. Seta Oxidation Vessels (15802-2) are recommended for use with the Seta Oxi Block Baths.

Pressure in the vessel can be monitored automatically with the Seta Auto-Oxi Control System.



## Specifications:

Bath Type:	Dry
Temperature Range:	Ambinet 0 to 130°C
Number of Test Stations:	4
Heat Up Time to 100°C	80 mins
Voltage:	200/240V 50/60Hz
Power Consumption:	3kW max
Size (HxWxD):	30x60x28 cm
Weight:	65 kg