

# Water Bath 90200-5

ASTM D96; ASTM D5546; API 2542; API 2548

15 litre capacity stainless steel water bath with Thermostir 120 unit mounted on top for pre-heating samples.

The Thermostir 120 unit comprises a 4 digit display digital temperature controller and stirrer in a single unit. The unit is designed to operate between ambient +5 °C to 120 °C with a setting resolution of 0.1 °C.

The water bath requires a stainless steel rack to hold the test tubes.

Supplied with: lid, mains lead, instruction manual.



## Specifications:

Temperature Range:	Ambient +5 to 120 °C
Temperature Stability:	±0.05 °C
Temperature Uniformity:	±0.01 °C
Display:	4 Digit LED
Setting Resolution:	0.01 °C
Size (HxWxD):	30 x 35 x 35 cm
Internal Bath Size (HxWxD):	33 x 30 x 15 cm
Weight:	8 kg