

# Seta Cloud and Pour Point Cryostats

ASTM D97; D2500; D5853; IP 15; 219; 441; BS 2000-15; 2000-219; BS EN 23015; ISO 3015; 3016; Pour Point of Petroleum Products; Cloud Point of Petroleum Products

- Three or four individually temperature controlled compartments
- Four air wells in each compartment
- CFC free refrigeration system
- Heated anti-condensation lid



## Ordering information

93531-7:	Seta Cloud and Pour Point Cryostat
94100-3:	Seta Compact Cloud and Pour Point Cryostat

Specifications:	93531-7	94100-3
Temperature range:	Ambient to -35 °C (x3) -35 to -51 °C (x1)	Ambient to -34 °C
Compartment volume:	2 Litres	2 Litres
Test positions:	16	12
Cool down time:	Approx 2 hrs (from 32 °C to -51 °C)	Approx 1 hr (from 32 °C to -34 °C)
Power supply:	1.8 kW 110/220 Vac, 50/60 Hz	750 W 220/240 Vac, 50 Hz
Size (HxWxD) / Weight:	950 x 630 x 620 mm / 115 kg	600 x 600 x 850 mm / 115 kg