

## Laboratory Refrigerator 93700-0

- Temperature range 0 to 10°C
- Cool storage of volatile materials
- CFC free
- Spark free

Laboratory refrigerator for safe cooling and chilled storage of volatile products, 'spark free' refrigerators are essential. The Seta laboratory refrigerator is of approved design with all electrical components isolated from the interior of the cabinent to reduce fire risk from internal vapours.

The storage compartment can be maintained at 0°C and is equipped with two storage shelves and a lower drawer. To safeguard the integrity of samples, the refrigerator is fitted with a lockable door.



93700-0 Specifications:	
Temperature Range:	0 to 10°C
Capacity:	150 litres
Voltage:	220/240V, 50/60Hz
Power:	70W
Size (HxWxD):	85x50x56cm (external) 75x45x45cm (compartment)
Weight:	20kg