

SetaCool 94150-0

Cloud, Pour Point and CFPP

ASTM D97; D5853; IP 15; 441; BS 2000-15; ISO 3016; Pour Point of Petroleum Products

ASTM D2500; IP 219; BS EN 23015; BS 2000-219; ISO 3015; Cloud Point of Petroleum Products

ASTM D6371; IP 309; EN ISO 116; Cold Filter Plugging Point of Diesel and Heating Fuels

- Built in cooling
- Dry block bath with no cooling medium required
- Independent modular system
- Low power consumption and energy loss
- Built-in heated top plate to prevent ice formation
- Temperature stability +/- 0.2 °C or better
- Small, compact footprint
- Low maintenance
- Eliminates greenhouse gas refrigerants



Technical Specifications	
Number of test jackets	4
Cooling temperature range	Adjustable between 0 to -69 °C with PID control
Temperature stability	+/- 0.2 °C or better
Power	100 V to 240 Vac 50/60 Hz
Size (HxWxD)	280 x 250 x 500 mm
Weight	20 kg