

Cold Filter Plugging Point Apparatus 99000-3

ASTM D6371; IP 309; BS EN 116; BS 2000-309

- Compatible with fuels containing flow improving additives

The Cold Filter Plugging Point (CFPP) Apparatus determines low temperature flow characteristics of diesel fuels, gas oils and other distillate fuels.

It comprises a 4 litre dry ice bath (Dewar Flask) and a CFPP test insert connected via a three way cock to a UKAS certified digital manometer and two 5 litre pressure equalising vessels.



Specifications:

Weight:	3.5 kg
---------	--------