

Items required for each test		ASTM D1657 (obs)	IP 59 (obs)	IP 235	ISO 3993
22580-2	Vapour Pressure Gauge 0 to 2000 kPa	●	●	●	●
22680-0	Hydrometer 0.50-0.58 gcm ³		●	●	●
22690-0	Hydrometer 0.57-0.65 gcm ³		●	●	●
22730-2	Thermometer: ISO653 -15 to 45 °C				●
22720-2	Thermometer: BS593 -20 to 60 °C			●	
22705-2	Hydrometer Thermometer ASTM310H	●			
22710-2	Hydrometer Thermometer ASTM101H	●			