


Seta Penetration Needles

18500-0 Seta Penetration Needle - Weight $2.5 \text{ g} \pm 0.05 \text{ g}$

Note: 18512-0 Seta Penetration Needle has the same dimensions and is UKAS certified to IP 49.




[ASTM: 50 to 55 mm (350 - 500 PEN)]

IP: > 40 to 45 mm and ferrule must not penetrate (330-500 x 0.1 mm)]

18520-0 Seta Standard Penetration Needle - Weight $2.5 \text{ g} \pm 0.05 \text{ g}$

Note: 18510-0 Seta Standard Penetration Needle has the same dimensions and is UKAS certified to IP 49.



[ASTM: 40 to 45 mm ($\leq 350 \times 0.1 \text{ mm}$)

IP: Approx. 40 to 45 mm ($\leq 330 \times 0.1 \text{ mm}$)]

18490-0 Seta Wax Penetration Needle - Weight $2.5 \text{ g} \pm 0.05 \text{ g}$

Note: 18492-0 Wax Penetration Needle has the same dimensions and is UKAS certified to ASTM D1321.

