

Items required for each test		Copper Corrosion		Silver Corrosion	
		ASTM D130	ASTM D4048	ASTM D4814	
		IP 154	IP 112	IP PM DJ 2005	IP 227
11241-0	Silicone Carbide Paper 150 grit (pack of 50)	●	●		●
11420-0	Support Assembly (pair)		●		
11430-0	Glassware Silver Corrosion				●
11440-0	Silver Test Strip (pack of 5)			■	●
11450-0	ASTM Silver Strip Corrosion Standard				●
11460-0	Silicon Carbide Paper, 65 micron, P240 FEPA Grade (pack of 50)	●	●	■	●
11470-0	Silicon Carbide Paper P100 FEPA Grade (pack of 50)	●	●		
11480-0	Silicon Carbide Grains 150 Mesh 500 g (pack of 5)	●	●	■	●
11490-0	Silver Strip Glass Cradle (pack of 3)			■	
11500-0	Pressure Vessel	●		■	
11510-0	Silver Corrosion Kit				▲
11514-0	Tongs	●	●	■	●
11520-0	Test Tube Ejector Spring (one per tube)			■	
11530-0	Multi-Strip Vice	●	●	■	●
11550-0	Copper Test Strip (pack of 30)	●	●		
11570-0	Flat Glass Test Tube (pack of 10)	●	●	■	
11580-0	ASTM Copper Strip Corrosion Standard	●	●	●	
11590-0	Test Tube (pack of 10)		●	■	
22750-3	Digital Stopwatch			●	
ASTMS12C	Thermometer	●	●	●	●

Items marked ■ are available as a kit (11515-0) for tests to ASTM D4814
 Item marked ▲ is only suitable for use with 15550-3 and 11400-7