

Items required for each test		Fuels	Oils	Oils	Oils
		*ASTM D2274	ASTM D943	ASTM D4310	DIN 51 587
11470-0	Silicon Carbide Paper P100 FEPA Grade (pack of 50)		●	●	●
16910-0	Set Glassware ASTM D2274, D7462, ISO 12205, IP 388 (6 cells)	●			●
16914-0	Sample Filtering Kit for ASTM D2274	●			
16915-0	Filtration & Accessory Kit (ASTM D4310)			●	
16920-0	Steel Catalyst Wire (pack of 6)		●	●	●
16921-0	Pre-Formed Copper and Steel Catalyst Coil (pack of 6)		●	●	●
16930-0	Copper Catalyst Wire (pack of 6)		●	●	●
16950-0	ASTM D943 & D4310 Accessory Pack		●		
16980-2	Thermometer for ASTM D943 & D4310 Oxidation Cell: Low Haz		●	●	
25000-4	Hot Plate	●			●
99100-2	Analytical Balance	●		●	●
99200-3	Laboratory Drying Oven	●		●	●

\* Note ASTM D2274 specifies the use of a liquid bath